



**CATALOGUE #:** 3NS8

**PRODUCT NAME:** Monoclonal mouse anti- Influenza A virus non structural (NS) protein

---

**MAbs:** 4A1, 9F10

**Specificity:** Non structural protein of influenza A virus.

**MAb isotypes:** **IgG2a** for MAb 4A1  
**IgG2b** for MAb 9F10

**Applications:** All MAbs can be used in ELISA.

**Presentation:** PBS, pH 7.4, 0.1 % sodium azide (NaN<sub>3</sub>)

**Storage:** + 4 °C

**Material safety note:** This product is sold as an antibody preparation **for research use only**. Standard Laboratory Practices should be followed when handling this material.

Contains sodium azide (0.1 %) as preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.