



CATALOGUE #: 4IG8

PRODUCT NAME: Monoclonal mouse anti-insulin-like growth factor binding protein-1 (IGFBP-1)

MAbs: G5F8, C7B9

Specificity: MAbs do not cross-react with human insulin or proinsulin

MAb isotypes: **IgG1** for MAb C7B9
IgG2a for MAb G5F8

Applications: Antibodies are working in ELISA and Western blotting.
Recommended pair for sandwich ELISA (coating – conjugate):
C7B9 – G5F8

Purification: Salt fractionation and ion-exchange chromatography

Presentation: PBS, pH 7.4, 0.04 % sodium azide (NaN₃)

Storage: + 4 °C

Material safety note: This product is sold as an antibody preparation **for research use only**. Standard Laboratory Practices should be followed when handling this material.

Contains sodium azide (0.04 %) as preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.