



CATALOGUE #: 4NN2

PRODUCT NAME: Monoclonal mouse anti-non-neuronal enolase (NNE)

MAbs: H11

Hybridoma clones have been derived from hybridization of Sp2/0 myeloma cells with spleen cells of Balb/c mice immunized with recombinant human NNE.

Specificity: NNE, no cross-reactivity with NSE

MAb isotypes: IgG1 for MAb H11

Applications: EIA, Western blotting

Purification: Chromatography on protein G Sepharose

Presentation: PBS, pH 7.4, 0.1 % sodium azide (NaN₃)

Storage: + 4 °C

Material safety note: This product is sold as an antibody preparation **for research use only**. Standard Laboratory Practices should be followed when handling this material.

Contains sodium azide (0.1 %) as preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.